

Purge&Trap Extraction of Volatiles in Drinking Water and High Speed TOF MS Detection in compliance to EPA 524.3

APPLICATION NOTE

Authors:

Ilaria Ferrante
Roberta Lariccia

DANI Instruments SpA
viale Brianza, 87
Cologno Monzese
Milano
Italy

Key Words:

EPA 524.3
Drinking Water
P&T
Fast GC-TOF MS

Introduction

Analysis of Volatiles Organic Compounds (VOCs) in water is of increasing interest due to critical ground contamination from several sources as gasoline, oil spills and industrial solvents.

Especially for drinking water, possible contamination is of great concern for potential human health effect. The U.S. Environmental Protection Agency (EPA) strictly regulates the assessment of drinking water quality through method 524.3 in which detection limits and instrumentation requirements are established. Official guidelines require the monitoring of VOCs contaminants in drinking water at progressively lower concentration level and typically the P&T extraction technique is indicated to reach the requested limit of detection.

This work is presenting the use of the flexible DANI Master Purge&Trap Sampler (Master P&T), highlighting the several benefits of this SYSTEM in

terms of extended automation, overlapping incubation time capability and absence of cross-contamination between samples.

Besides, particular attention has been paid to the most volatiles compounds in terms of peak shape and recovery.



Figure 1 : Master P&T coupled to the Master GC-TOF MS

Environmental

DANI

A SCENT OF FUTURE

Method

In this work the Purge&Trap extraction method, coupled with a GC-TOF MS analysis, is described to determine VOCs in water.

The system is illustrated in Figure 1. Coupling the Purge&Trap with an automatic sampler guarantees a good repeatability in terms of internal standard addition, assuring a complete automation of all operation steps.

In Figure 2 the schematic plumbing of the P&T is illustrated. The use of a triple layer sorbent trap has been found to be ideal to handle the wide range of volatility and polarity represented by the compounds under study.

The use of the Master TOF MS as detector offers the possibility to discriminate even between coelute compounds exploiting the Deconvolution Algorithm of MasterLab Processing software.




Figure 2: Dynamic Headspace and Purge&Trap Sampler schematic plumbing in the "trapping" phase

Sample Preparation

Calibration standard mixtures and fortified blanks have been prepared following the method guideline, starting from commercial certified VOCs mixture in methanol as stock solution.

A Primary Dilution Standard in methanol has been used to spike 10 mL of organic-free water (HPLC grade) directly into the vial to prepare calibration standards for linearity evaluation in the range of 0.02 to 50 µg/L. Internal Standard required from EPA method 524.3 were added in a final concentration of 2 ppb to each sample, using MasterAS.

According to quality control criteria, before

calibration, the TOF-MS has been checked to meet tuning requirement by injecting 4-bromofluorobenzene (BFB). All the ion ratio criteria is achieved after a standard tuning of the system. (Figure 8)

| Master GC | |
|-----------|---|
| Column | wozol 20 m, 0.18 mm, 1µm |
| Flow Rate | 0.8 mL/min |
| SSL Temp | 250°C |
| Oven Temp | 35 (5 min) to 240°C (2 min) at 20°C/min |

| Master TOF | |
|--------------------|-------------|
| Transfer Line Temp | 250°C |
| EI source | 70 eV |
| Source Temperature | 300°C |
| Mass Range | 35-360 amu |
| Acquisition Rate | 5 spectra/s |

| Master DHS/P&T in Purging Mode | |
|--------------------------------|------------------------------|
| Sample Volume | 5 mL |
| Incubation Temp | 60°C |
| Stripping | 5 min @ 100 mL/min at Helium |
| Trap Temp | -10°C |
| Oven Stop Temp | 0°C |
| Trap Desorption Temp | 200°C |
| Injection Time | 2 min |
| Transfer Line Temp | 250°C |
| Valve Temp | 250°C |
| Baking Time | 10 min @ 80 mL/min |
| Trap Baking Temp | 205°C |
| Oven Stop Baking Temp | 100°C |
| Trap Sorbent | Tenax/ Carbotrap/ Carboxene |

Table 1: Operative Conditions

Results and Discussion

Faster separation has been obtained by using a shorter and narrow bore column compared to conventional approach using longer column. In Figure 4 the TIC of a mixture containing 52 compounds (listed in Table 2) is reported. Using a 20m column with an id of 0.18mm it is possible to complete the separation in about **14 minutes**, however several compounds appear partially or completely coeluted. The advantage of using the TOF MS detector is the possibility to discriminate within coeluted compounds due to powerful deconvolution capability. With DANI MasterLAB processing software (MLP) it is possible to locate

different apexes under the same TIC peak, allowing the proper identification of coeluted compounds. As an example, Figure 4 shows the deconvolution results for the peak at RT 9.93 min, with the identification of two different compounds.

| RT | Compound Name | Integration Source | Linearity r ² | MDL (ppt) | EPA Limits (ppt) |
|--------|--------------------------------|----------------------|--------------------------|-----------|------------------|
| 4.528 | Glyoxal, 1,1-dichloro- | m/z 61 (See Note 6) | 0.998 | 9.94 | 40 |
| 4.630 | Glyoxal, 1,1-dichloro- | m/z 62 (See Note 6) | 0.993 | 1.91 | 20 |
| 7.017 | cis-2,3-dichlorocyclohexane | m/z 61 (See Note 6) | 0.999 | 2.40 | 42 |
| 7.288 | Methane, bromochloro- | m/z 150 (See Note 6) | 0.999 | 2.67 | 42 |
| 7.327 | Trichloromethane | m/z 62 (See Note 6) | 0.999 | 0.50 | 25 |
| 7.528 | Glyoxal, 1,1,1-trichloro- | m/z 97 (See Note 6) | 0.998 | 1.21 | 26 |
| 7.687 | 1-Propene, 1,1-dichloro- | m/z 75 (See Note 6) | 0.999 | 2.15 | 82 |
| 7.725 | Carbon tetrachloride | m/z 117 (See Note 6) | 0.998 | 2.88 | 44 |
| 7.988 | Benzene* | Ion (deconvolution) | 0.999 | 6.41 | 17 |
| 7.912 | Glyoxal, 1,2-dichloro* | Ion (deconvolution) | 0.993 | 7.54 | 25 |
| 8.218 | 1,2-Benzenediol, 4,4-dichloro- | m/z 114 (See Note 6) | IS | IS | IS |
| 8.655 | Trichloroethylene | m/z 129 (See Note 6) | 0.999 | 1.11 | 25 |
| 8.632 | Propene, 1,2-dichloro- | m/z 62 (See Note 6) | 0.999 | 1.88 | 18 |
| 8.792 | Methane, dibromo- | m/z 174 (See Note 6) | 0.999 | 2.45 | 45 |
| 8.888 | Methane, bromodichloro- | m/z 62 (See Note 6) | 0.999 | 0.95 | 14 |
| 9.228 | cis-2,3-dichloropropene | m/z 75 (See Note 6) | 0.999 | 2.16 | 26 |
| 9.488 | Toluene | m/z 91 (See Note 6) | 0.999 | 0.25 | 24 |
| 9.642 | trans-1,2-dichloropropene | m/z 75 (See Note 6) | 0.999 | 0.94 | 12 |
| 9.795 | Glyoxal, 1,1,2-trichloro- | m/z 97 (See Note 6) | 0.987 | 0.89 | 48 |
| 9.915 | Propene, 1,2-dichloro* | Ion (deconvolution) | 0.999 | 5.83 | 20 |
| 9.928 | Trichloroethylene* | Ion (deconvolution) | 0.998 | 1.28 | 26 |
| 10.118 | Methane, dibromodichloro- | m/z 129 (See Note 6) | 0.999 | 5.87 | 27 |
| 10.288 | Glyoxal, 1,2-dibromo- | m/z 187 (See Note 6) | 0.998 | 2.98 | 18 |
| 10.542 | 1,2-Dichlorobenzene-d5 | m/z 117 (See Note 6) | IS | IS | IS |
| 10.562 | Benzene, chloro- | m/z 112 (See Note 6) | 0.998 | 0.58 | 19 |
| 10.628 | Glyoxal, 1,2-dichloro- | m/z 91 (See Note 6) | 0.999 | 0.56 | 10 |
| 10.712 | m-ly-Styrene | m/z 91 (See Note 6) | 0.993 | 1.28 | 20 |
| 11.012 | o-Styrene* | Ion (deconvolution) | 0.998 | 6.88 | 10 |
| 11.027 | Styrene* | Ion (deconvolution) | 0.998 | 4.48 | 11 |
| 11.278 | Isopropylbenzene | m/z 91 (See Note 6) | 0.987 | 0.87 | 11 |
| 11.580 | Glyoxal, 1,1,2,2-tetrachloro- | m/z 131 (See Note 6) | 0.998 | 5.41 | 13 |
| 11.528 | Benzene, bromo* | Ion (deconvolution) | 0.998 | 28.15 | 20 |
| 11.547 | Propene, 1,2,3-trichloro* | Ion (deconvolution) | 0.998 | 17.28 | 50 |
| 11.582 | n-propylbenzene | m/z 91 (See Note 6) | 0.998 | 0.76 | 77 |
| 11.688 | 2-chlorocyclohexane | m/z 91 (See Note 6) | 0.998 | 4.46 | 21 |
| 11.797 | Benzene, 1,2,4-trichloro- | m/z 135 (See Note 6) | 0.987 | 0.88 | 15 |
| 11.792 | 4-chlorocyclohexane | m/z 91 (See Note 6) | 0.987 | 0.85 | 14 |
| 11.885 | Benzene, iso-butyl- | m/z 119 (See Note 6) | 0.998 | 0.54 | 20 |
| 11.888 | Benzene, 1,2,4-trichloro- | m/z 135 (See Note 6) | 0.998 | 0.55 | 80 |
| 12.117 | sec-butylbenzene | m/z 135 (See Note 6) | 0.998 | 1.71 | 12 |
| 12.237 | 4-isopropyltoluene | m/z 119 (See Note 6) | 0.998 | 0.71 | 12 |
| 12.227 | Benzene, 1,4-dichloro- | m/z 146 (See Note 6) | 0.998 | 2.62 | 15 |
| 12.270 | 1,2,4-Trichlorobenzene-3H | m/z 159 (See Note 6) | IS | IS | IS |
| 12.288 | Benzene, 1,3-dichloro- | m/z 146 (See Note 6) | 0.998 | 2.79 | 12 |
| 12.525 | Benzene, butyl- | m/z 91 (See Note 6) | 0.998 | 1.45 | 45 |
| 12.577 | Benzene, 1,2-dichloro- | m/z 146 (See Note 6) | 0.998 | 1.34 | 19 |
| 13.185 | Propene, 1,2-dibromo-1-chloro- | m/z 157 (See Note 6) | 0.998 | 15.88 | 41 |
| 13.788 | Benzene, 1,2,4-trichloro- | m/z 139 (See Note 6) | 0.987 | 12.86 | 20 |
| 13.885 | 1,1,1,2,2,2-hexachloroethane | m/z 225 (See Note 6) | 0.998 | 12.18 | 42 |
| 13.875 | Nyloleone | m/z 129 (See Note 6) | 0.993 | 4.88 | 12 |
| 14.148 | Benzene, 1,2,4-trichloro- | m/z 139 (See Note 6) | 0.998 | 19.81 | 11 |

Table 2: List of analyzed VOCs and extrapolated MDLs (3x noise level)

Some of the compounds that EPA requires to keep under control are very volatile even at low temperature, such as the vinyl chloride that has a boiling point of -13.4°C. A solution containing only the very volatile compounds has been used to evaluate the repeatability of the injection (Table 3). The linearity was very good over the required concentration range with correlation coefficients > 0.995 for all the compounds as listed in Table 2.

| Compound Name | Base Peak | Average Area | SD | RSD% |
|--------------------------|-----------|--------------|---------|------|
| Methane, chloro- | m/z 50 | 13283.67 | 300.59 | 2.26 |
| Ethane, chloro- | m/z 62 | 34586.67 | 2001.15 | 5.79 |
| Methane, bromo- | m/z 94 | 12480.67 | 389.91 | 3.04 |
| Ethyl Chloride | m/z 64 | 12825.67 | 389.91 | 3.04 |
| Trichloromethane methane | m/z 101 | 30880.00 | 1156.15 | 3.74 |

Table 3: Repeatability on 3 consecutive extraction of the most volatile compounds




Figure 3: Example of Linearity response




Figure 4 : Fast GC Chromatogram (TIC) of a 52 VOCs mixture and example of deconvolution elaboration (a)




Figure 5 : Fast GC Chromatogram (TIC) of a 52 VOCs mixture and example of deconvolution elaboration (a)




Figure 6 : Some compounds required by the EPA norm are very volatile even at low temperature. Among these, there is the Vinyl Chloride. Even for such compound, a good quantification has allowed to obtain a very good peak shape.





Figure 7 : Thanks to the good peak shapes, it is possible to easily calculate the calibration plot and obtain a very good minimum detectable level like 3 times the noise on the lower concentration of calibration (2 ppt)



| Ion m/z | EPA 524 Ion Ratio Criteria | MasterTOF % Relative Intensity | Pass/Fail |
|------------|-------------------------------|-----------------------------------|-----------|
| 50 | 15-40 % of mass 95 | 15.7 | ✓ PASS |
| 75 | 30-60 % of mass 95 | 39.1 | ✓ PASS |
| 95 | Base peak | 100 | ✓ PASS |
| 96 | 5-9 % of mass 95 | 7.1 | ✓ PASS |
| 173 | < 2 % of mass 174 | 1.8 | ✓ PASS |
| 174 | > 50 % of mass 95 | 74.6 | ✓ PASS |
| 175 | 5-9 % of mass 174 | 7.8 | ✓ PASS |
| 176 | 95-105 % of mass 174 | 97.0 | ✓ PASS |
| 177 | 5-10 % of mass 176 | 6.8 | ✓ PASS |

Figure 8 : EPA method 524.3 requires to check the quality of fragmentation for the BFB. All the ion ratio criteria can be achieved with a standard tuning of the system without any specific tuning required.



A SCENT OF FUTURE



A SCENT OF FUTURE

www.danipa.it

The contents of this publication are for reference and illustrative purposes only. Information, descriptions and specification are subjected to change without notice. DANI Instruments assumes no responsibility and will not be liable for any errors or omissions contained herein or for incidental, consequential damages or losses in connection with the furnishing, performance, or use of this material.

AN_131
