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

The average value of coconut juice was -25.4 $\delta^{13}\text{C}$, ‰, and the coconut water with sugar was -19.1 $\delta^{13}\text{C}$, ‰.

4CUS Service Info C 0 1 1

CAVITY RING - DOWN SPECTROSCOPY (CRDS)

CM-CRDS is a powerful analytical platform in the food and beverage adulteration market. It can easily distinguish down to ~5% adulteration of coconut water by the addition of C4 sugars, with excellent precision (about ± 0.3 per mil) and accuracy.

What is CRDS?

Nearly every small gas-phase molecule (e.g., CO_2 , H_2O , H_2S , NH_3) has a unique near-infrared absorption spectrum. At sub-atmospheric pressure, this consists of a series of narrow, well-resolved, sharp lines, each at a characteristic wavelength. CRDS - Cavity Ring-Down Spectroscopy - avoids this sensitivity limitation by using an effective pathlength of many kilometers. It enables gases to be monitored in seconds or less at the parts per billion level, and some gases at the parts per trillion level. The beam from a single-frequency laser diode enters a cavity defined by two or more high reflectivity mirrors. The analyzers use a three-mirror cavity to support a continuous traveling light wave. This provides superior signal to noise compared to a two-mirror cavity that supports a standing wave. When the laser is on,

In this note, we present the Combustion Module Cavity Ring-Down Spectroscopy (CM-CRDS) System as a low-cost, easy-to-use, accurate, and precise solution to identify the addition of C4 sugars (Adulteration), based on the carbon isotopic ($\delta^{13}\text{C}$) composition of the coconut water. It is easier to use and a significantly lower capital and total-ownership cost alternative to Isotope Ratio Mass Spectrometry (IRMS), and it produces results that are indistinguishable. The CM-CRDS system is simple to install, easy to operate, and does not require sample preparation.